

# Bacterial Transformation

(Chemically Competent Cells)

Leslie Vosshall

1/22/2001

- 1) Retrieve the required number of tubes of competent cells and place on wet ice. Our lab stock is XL1-Blue MRF<sup>7</sup> and is packaged in 400 ul aliquots stored at  $-80^{\circ}\text{C}$ . You will need 100 ul of competent cells for each transformation. Once thawed, aliquots cannot be re-frozen, so discard extra cell volume if you are only transforming 1 or 2 samples.
- 2) Thaw the competent cells on ice. You can accelerate the thawing by briefly warming the tubes in your hands.
- 3) Once competent cells are thawed, dispense 100 ul of competent cells into each sample, consisting of either 5 ul of a ligation or 0.5 ul of a pure plasmid to be retransformed. Flick tube gently to mix DNA + Competent Cells.
- 4) Incubate DNA + Competent Cells on wet ice for 20'-30'.
- 5) Heat shock at  $42^{\circ}\text{C}$  for EXACTLY 45". Do not exceed 45".
- 6) Allow cells to recover on wet ice for 2 minutes.
- 7) Add 200 ul of SOC to cells and incubate at  $37^{\circ}\text{C}$  for 20'-30'.
- 8) Briefly spin cells down in capsulefuge, discard all but ~100 ul of supernatant, resuspend cells and plate entire volume on LB plates including appropriate antibiotic. Invert plates and incubate overnight at  $37^{\circ}\text{C}$ .