

# Enhanced Two-Color Fluorescent *In Situ* Hybridization for Low-Abundance Genes

## ***In situ* hybridization two-color fluorescent detection using TSA-FITC/ TSA-CY5**

Leslie Vosshall 2 September 2000

Elane Fishilevich 21 April 2005

Revised: Kirsten Vannice 25 January 2006

### **DAY ONE**

\*Prepare fresh frozen tissue in OCT and section at 14  $\mu$ m. Collect sections on Superfrost PLUS slides. Let slides dry 20 minutes to 3 hours. \*Fix in 4% Paraformaldehyde/1XPBS 10 minutes

To make fixative:

Add 20g paraformaldehyde to 450 ml boiling dH<sub>2</sub>O + 50  $\mu$ l 10N NaOH. Stir on magnetic stirrer in hood until paraformaldehyde is well dissolved. Add 50 ml 10XPBS. Sterile filter and chill on ice. Use within 3 days of preparation.

\*Wash 3X3 minutes in 1XPBS.

\*Acetylate for 10 minutes.

To make acetylation solution:

Combine 9.25 g triethanolamine HCl Sigma Cat.#T1502 and 1.12 ml 10N NaOH in 500 ml dH<sub>2</sub>O. Mix well. Immediately before adding solution to sections, add 1.25 ml acetic anhydride and shake solution vigorously.

\*Wash 3X5 minutes in 1XPBS.

\*Prehybridize by placing slides vertically in slide mailers filled with hybridization solution at 65°C for 1 hour. You can reuse hybridization solution for 1-2 months if you store it at -20°C.

#### **Hybridization solution:**

50% formamide

5X SSC

5X Denhardt's

250  $\mu$ g/ml yeast tRNA

500  $\mu$ g/ml herring sperm DNA

50  $\mu$ g/ml Heparin

2.5 mM EDTA

0.1% Tween-20

0.25% CHAPS

\*For each slide, prepare 100 µl hybridization buffer + 1-2 µl of DIG and/or FITC probe (see section below for probe synthesis; can adjust amount of probe used based on strength, eg between 1:2500 and 1:10).

\*Denature probes by heating to 80 degrees C for 5 minutes. Chill on ice. Pipet carefully onto tissue and coverslip with Hybrislip plastic cover slips.

\*Place slides horizontally in slide mailers. Do not mix slides containing different probes in the same slide mailer, because cross-contamination may occur.

\*Hybridize slides by placing horizontally in slide mailers in a humidified box at 65oC for 16 to 24 hours.

## **DAY TWO**

Necessary Reagents:

anti-FLUORESCEIN-POD, Fab fragments (Roche 1 426 346)

anti-DIGOXIGENIN-POD, Fab fragments (Roche 11 207 733 910)

TSA™ Plus Fluorescein Fluorescence System (Perkin Elmer NEL741)

TSA™ Plus Cy5 Fluorescence System (Perkin Elmer NEL745)

Roche

Anti-digoxigenin-POD, Fab fragments

Catalog 11 207 733 910

Roche

Anti-Fluorescein-POD, Fab fragments

1 426 346

Buffers:

TN

100 mM Tris-HCl (pH 7.5)

150 mM NaCl

TNB

0.5% Blocking Reagent (Perkin Elmer) in TN buffer

TNT

0.05% Tween20® in TN buffer

\*Remove coverslips by soaking slides in 5X SSC prewarmed to 65oC

\*Wash 3X20 minutes in 0.2X SSC at 65oC

\*Equilibrate in TN for 5 min.

\*Block in TNB for 30 min in a vertical slide mailer @ RT.

\*Dilute anti-FITC-POD 1:500 in TNB.

Apply 100 µl/ slide, coverslip, and incubate for 30 min @ RT.

\*Wash slides with TNT - 3 x 5 min.

\*Dilute FITC-Tyramide 1:50 in Amplification Reagent (TSA™ kit).

Apply 100 µl of working solution per slide, coverslip, incubate for 10 min @ RT.

\*Wash slides with TNT - 1 x 5 min.

\*To quench peroxidase activity, incubate in 3% H<sub>2</sub>O<sub>2</sub> for 1hr.

\*Wash slides with TNT - 3 x 5 min.

\*Dilute anti-DIG-POD 1:500 in TNB.

Apply 100 µl/ slide, coverslip, and incubate for 30 min @ RT.

\*Wash slides with TNT - 3 x 5 min.

\*Dilute Cy5-Tyramide 1:50 in Amplification Reagent (TSA™ kit).

Apply 100 µl of working solution per slide, coverslip, incubate for 10 min @ RT.

\*Wash slides with TNT - 3 x 5 min.

\*Mount Slides in 60ul VectaShield®.

NOTE: If doing double in situ, do weaker probe first. Eg, if DIG probe is weaker than FITC probe, incubate with anti-DIG-POD first and anti-FITC-POD second. If using this protocol for a single in situ, ignore quenching.

## **DIG RNA RIBOPROBES**

Leslie Vosshall 2 September 2000

Revised: Kirsten Vannice 25 January 2006

In my experience, the purity of the template DNA is the major determinant of successful RNA synthesis. This protocol therefore takes special care in post-digestion purification. It is recommended that these steps are followed to ensure consistently good RNA probes.

### Make Template

\*Digest the equivalent of 1  $\mu\text{g}$  of insert DNA in 50  $\mu\text{l}$  with 3  $\mu\text{l}$  of enzyme (to generate blunt ends or 5' overhang; no 3' overhangs) for 3 hours. [Calculate amount of DNA to yield 1  $\mu\text{g}$  of insert. e.g. 0.5 kb insert in pBluescript--digest 7  $\mu\text{g}$  of plasmid DNA; **also note** .5 $\mu\text{g}$  of insert DNA can be used with success if DNA is limited]

\*Check 1  $\mu\text{l}$  of digest in a gel to make sure DNA is linearized.

\*If linear, dilute remaining 49  $\mu\text{l}$  of digest with 151  $\mu\text{l}$  RNase free water and purify DNA by extracting once with 200  $\mu\text{l}$  of 1:1 phenol:chloroform and once with 200  $\mu\text{l}$  of chloroform. Vortex each extraction for 3 minutes and spin for 4 minutes at maximum speed at room temperature.

\*Pass aqueous phase over a G50 ProbeQuant column (Pharmacia) to remove residual organic solvents.

\*Measure volume and precipitate eluate with 10%  $\mu\text{l}$  sodium acetate (3M, pH5) and 250%  $\mu\text{l}$  ethanol. Spin out DNA, 30 minutes 4 degrees. Wash pellet with 70% ethanol. Dry and resuspend in 11  $\mu\text{l}$  of RNase free water

\*Check .5  $\mu\text{l}$  of DNA on gel to make sure template is still present.

### Synthesize Riboprobes

\*To 10.5  $\mu\text{l}$  of template in dH<sub>2</sub>O

add in order for a total volume of 20  $\mu$ l:

4  $\mu$ l Promega 5X Transcription Buffer

2  $\mu$ l BMB 10X DIG RNA labeling mix (BMB Cat. #1 277 073)

1  $\mu$ l 0.1M DTT

0.5  $\mu$ l Promega rRNAsin

2  $\mu$ l RNA polymerase (SP6, T3, or T7)

\*Mix gently, briefly spin down tube at low speed, and incubate at 37 degrees for 2 hours.

\*Pass probes over G50 ProbeQuant columns (Amersham/Pharmacia Cat.#)

\*Check 2  $\mu$ l of probe on a gel (can use normal agarose gel) to ensure transcription reaction was successful.

\*Bring volume of eluate to a final volume of 100  $\mu$ l with hybridization buffer (i.e., add ~82  $\mu$ l). Use 1 $\mu$ l-2 $\mu$ l/ml for in situ hybridization.

#### FITC RNA RIBOPROBES

Substitute 10X Fluorescein RNA labeling mix (BMB Cat. #1 685 619) for 10X DIG RNA labeling mix in protocol above. Do all incubations in the dark. Eluted probe will be yellow if synthesis has succeeded.

#### BIOTIN RNA RIBOPROBES

Substitute 10X BIOTIN RNA labeling mix (BMB Cat. #1 685 597) for 10X DIG RNA labeling mix in protocol above.