

Genomic Southern Protocol

6/27/2002 Leslie Vosshall

- 1) Cut 10 fly equivalents of DNA (or 10 ug of mouse/rat/human genomic DNA) with 5 ul appropriate enzyme in 100 microliter total volume, overnight.
- 2) Ethanol precipitate digest (1/10 volume 3M NaOAc pH 5.2, 2.5 volumes absolute ethanol), -20 degrees C for 30 minutes. Spin down DNA, 30 minutes, 15,000 rpm, 4 degrees C. Carefully remove supernatant, add 1 ml 70% ethanol. Vortex 10 seconds. Spin 5 minutes, 15,000 rpm. Carefully remove supernatant, dry down pellet in speedvac briefly. Do not overdry!
- 3) Resuspend DNA in 30 or 40 microliters TE (10mM Tris HCl pH 7.5-8.0, 1mM EDTA), depending on whether you will run a 15 or 20 well gel. If using Owl 15 well comb=55ul/well final, If using Owl 20 well comb=40ul/well final volume. Make sure that DNA is completely resuspended by gently pipeting up and down and heating sample at 55 degrees C and then gently mixing. The pellet from ethanol precipitation should be completely dissolved.
- 4) Pour a 0.8% agarose gel, casting it with a single 15 or 20 well comb. Plan on having one marker lane on each end of the gel with your samples in between. Pour a sufficient number of gels to accommodate all your samples. Depending on size of bands expected, load either 1KB Plus ladder and/or lambda/HindIII digested DNA as molecular weight marker. Add 5 ul loading dye to sample and heat, 65 degrees C for 5 minutes. Load gel and run slowly, 50 V for 6 hours or 20 mAmps overnight.
- 5) Take a picture of gel with fluorescent ruler on either side, aligned with zero at wells, so you can align markers later. Let gel linger under UV illumination on transilluminator for a few minutes to nick DNA (helps transfer of high molecular weight DNAs). Trim gel by cutting off the lower left corner, trim off wells.
- 6) Acid treat gel 15 minutes at room temperature with 12.5 ml glacial acetic acid/liter of dH2O with agitation.
- 6) Rinse gel with dH2O and incubate with denaturing solution [0.5N NaOH/1.5M NaCl] with agitation for 30 minutes.
- 7) Remove denaturing solution, rinse with dH2O, incubate gel in neutralization solution with agitation for 30 minutes to 60 minutes [1.5M NaCl/1M Tris pH 7.5].
- 8) Remove neutralization solution, incubate gel in 20X SSC for 5 minutes.
- 9) Set up downward transfer apparatus (Turbo Blotter from Schleicher and Schuell), using 20X SSC as transfer buffer and Hybond N+ as transfer membrane. If you plan to do low stringency hybridization with formamide, use Nitropure (MSI) not Hybond N+ (Amersham). Transfer for 3 hours only!
- 10) Disassemble transfer and rinse filter with 2XSSC. Mark wells with a black ballpoint pen. UV crosslink DNA to membrane with UV crosslinker. If using Nitropure bake in vacuum oven for 2 hours. Keep filter in 2XSSC until ready to hybridize, at 4 degrees if delay is longer than a few hours.
- 11) Prehybridize in Buck Buffer (20-40 ml) 1 hour to 6 hours, 65 to 70 degrees C. (0.5M Sodium Phosphate pH 7.3, 1% BSA, 4% SDS, heated to 65 degrees C if precipitates are present. Add sheared salmon sperm DNA, boiled for 10 minutes, then ice chilled to a final concentration of 0.5 mg/ml)
- 12) Remove prehybridization solution, denature probe, and add to more Buck Buffer (10-20 ml). Hybridize overnight at 65 to 70 degrees C.(~16 hours).
- 13) Remove probe and freeze -20 degrees C. Can reuse within one week.
- 14) Wash 3 X 30 minutes with prewarmed 0.2XSSC/0.5% SDS, 65 to 70 degrees C., with agitation.
- 15) Remove membrane from wash solution, blot dry on filter paper, tape filter to 3mm paper, orient with fluorescent ink, and cover with saran wrap. Expose to XOMAT-AR X-Ray film in cassette with intensifying screen, -80 degrees C, overnight. Develop film and put up for a longer exposure as necessary.